



Approved Cleaners August 2019

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

MANUFACTURER	PRODUCT	MAIN ACTIVE COMPONENT	CONCENTRATION	DATE TESTED
Ecolab®	A-456-II		0.4% A-256-II (1:256 Water)	02-Dec-14
Ecolab®	Actichlor Plus	Sodium Dichloroisochyanurate 30-60%	0.2% Actichlor Plus (2000PPM)	20-Nov-12
Quantum Pacific Limited	Actimax		According to Manufacturer Spec	11-Apr-06
SCIENCECompany	Ammonium Chloride	Cocoalkyldimethylbenzyl ammonium chloride	According to Manufacturer Spec	17-May-19
Milton®	Anti-Bacterial Tablets	Sodium Dichloroisochyanurate 30-60%	1 Tablet:2L Water	19-Jul-12
Ecolab®	Asepti Active Liquid		According to Manufacturer Spec	20-Nov-12
Ecolab®	Asepti Active Wipes	Dual-Quaternary/Alcohol	Ready to Use	20-Nov-12
Synergy Health	Azomax™ Multi Surface Wipes	Benzylkoniumchloride 21.6%	Ready to Use	15-Mar-13
Synergy Health	Azowipe™ Hard Surface Wipes	Benzylkoniumchloride 2%	Ready to Use	15-Mar-13
PreVasive®	BAC Botanical Antimicrobial Cleaner	Thyme Oil	According to Manufacturer Spec	22-May-19
Spartan Chemical Company Inc.	BNC-15™		According to Manufacturer Spec	11-Aug-16
Metrex® Research	CaviCide®	Isopropanol	According to Manufacturer Spec	06-Oct-16
Metrex® Research	CaviCide1™	Isopropanol	According to Manufacturer Spec	05-Jan-17
Metrex® Research	CaviWipes1™	Isopropanol	Ready to Use	05-Jan-17
Metrex® Research	CaviWipes™	Isopropanol	Ready to Use	06-Oct-16



Approved Cleaners August 2019

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

MANUFACTURER	PRODUCT	MAIN ACTIVE COMPONENT	CONCENTRATION	DATE TESTED
Agar Cleaning Systems Pty. Ltd.	Chloradet	Sodium Hypochlorite <10%	0.1% Cloradet (1:1000 Water)	06-Jun-16
Apollo Corporation	CID A-L™ II Quaternary Disinfectant	Quaternary Ammonium Compounds	According to Manufacturer Spec	26-Sep-12
Wexford Labs, Inc.	CleanCide® Wipes	Citric Acid	Ready to Use	15-May-19
Clorox Disinfecting Wipes	Clorox Disinfecting Wipes		According to Manufacturer Spec	14-Feb-11
Alco Mandala Products	DC-7		According to Manufacturer Spec	16-Feb-12
Clinell®	Detergent Wipes	Propane 1,2-diol	According to Manufacturer Spec	27-May-09
ArjoHuntleigh	Disinfectant Cleanser IV		According to Manufacturer Spec	23-Nov-09
Clorox®	Dispatch® Disinfectant		According to Manufacturer Spec	16-Feb-12
Tristel Solutions LTD	Distel Surface Disinfectant (Formerly Trigene®)	Polymeric Biguanide Hydrochloride <10%	According to Manufacturer Spec	24-Oct-06
Sealed Air Diversey Care	Divercleanse	Sodium Hypochlorite <5%, Sodium Peroxide <2%	3.1% Divercleanse (3.1:100 Water)	15-Dec-09
Buckeye® International, Inc.	Eco E23 S23		256 to 1 dilution with water	05-Jan-17
Velocity Chemicals	Eco-Pine		5% Ecopine (1:20 Water)	11-Apr-06
Clorox HEALTHCARE®	Fuzion™		According to Manufacturer Spec	05-Jan-17
Spartan Chemical Company Inc.	GS Neutral Disinfectant #3502		According to Manufacturer Spec	02-Dec-14
Clinell®	Hand & Surface Sanitizer	Benzalkonium Chloride	According to Manufacturer Spec	27-May-09



Approved Cleaners August 2019

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

MANUFACTURER	PRODUCT	MAIN ACTIVE COMPONENT	CONCENTRATION	DATE TESTED
Spartan Chemical Company Inc.	hdq®C2	Ammonium Chloride	According to Manufacturer Spec	30-Jun-17
Wipesplus	Hospital/Surgical Wipes	Ammonium Chloride	Ready to Use	15-May-19
Clorox HEALTHCARE®	Hydrogen Peroxide Cleaner Disinfectant Wipes	Hydrogen Peroxide	Ready to Use	12-Apr-19
Apollo Corporation	Hygena™	Ethylalcohol 1.4%	According to Manufacturer Spec	25-Sep-12
Kimberly-Clark®	Isowipes®	Isopropanol 30-60%	Ready to Use	19-Jul-12
Legge Systems	Legcide HB Plus		According to Manufacturer Spec	*****
Whiteley Medical	Lemex		According to Manufacturer Spec	15-Mar-11
Livingstone International (PTY) LTD	LIV-Wipes Jumbo	70% Isopropyl Alcohol	Ready to Use	15-Mar-11
Sulco LTD	Mediwipes		Ready to Use	27-May-09
Medline Industries Inc.	Micro-Kill™ Q3 Concentrated Disinfectant, Cleaner & Deodorizer		According to Manufacturer Spec	01-Feb-18
Medline Industries Inc.	MICRO-KILL™ AF2 Wipes	Ammonium Chloride	Ready to Use	15-May-19
Medline Industries Inc.	MICRO-KILL™ Bleach-Free, Alcohol-Free Wipes	Ammonium Chloride	Ready to Use	15-May-19
Medline Industries Inc.	MICRO-KILL™ One Wipes	Ethyl Alcohol	Ready to Use	15-May-19
Pro Restore® Products	Microban X-590 Institutional Spray	Diisobutylphenoxyethoxy ethyl dimethyl benzyl ammonium chloride	According to Manufacturer Spec	11-Aug-16
Anti-Slip Anywhere	Multi Kleen		According to Manufacturer Spec	15-Mar-11



Approved Cleaners August 2019

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

MANUFACTURER	PRODUCT	MAIN ACTIVE COMPONENT	CONCENTRATION	DATE TESTED
Certo®	Non Acid Bathroom Cleaner		According to Manufacturer Spec	07-May-12
PreVasive®	Noroxycdiff™	Hydrogen Peroxide	According to Manufacturer Spec	22-May-19
Ecolab®	Oasis Pro 12 Neutral Cleaner	Alcohols <10%	5% Oasis Pro 12 (1:20 Water)	06-Mar-12
Whiteley Medical	Odacon Spray		According to Manufacturer Spec	14-Jun-11
OdorXit Naturals	Odor Eliminator Spray		Ready to Use	22-Feb-18
Pal® International	PAL® Detergent Wipes		Ready to Use	27-May-09
Pal® International	PAL® Disinfectant Wipes		Ready to Use	27-May-09
PCS Process Cleaning Solutions	PCS 7000	Sodium Hypochlorite	According to Manufacturer Spec	06-Oct-17
PCS Process Cleaning Solutions	PCS Spray & Wipe Hand and Surface Cleaner		Ready to Use	06-Oct-17
Pure Green, LLC	Puregreen 24 (Pure Green)		According to Manufacturer Spec	30-Jan-09
Reynard Health Supplies	Reynard Neutral Detergent Wipes		According to Manufacturer Spec	08-Feb-07
Reynard Health Supplies	Reynard Surface Sanitizing Wipes		According to Manufacturer Spec	27-May-11
Sealed Air Diversey Care	Room Care R2L		According to Manufacturer Spec	24-Oct-12
STERI-7®	S-7 Xtra Wipes		Ready to Use	28-Nov-16
Alco Mandala Products	Sanni-Rinse		According to Manufacturer Spec	16-Feb-12



Approved Cleaners August 2019

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

MANUFACTURER	PRODUCT	MAIN ACTIVE COMPONENT	CONCENTRATION	DATE TESTED
Clinell®	Sporicidal Wipes		According to Manufacturer Spec	27-May-09
Ecolab®	Spring Clean		5% Spring Clean (1:20 Water)	11-Apr-06
Noble Pine Products Co.	Steri-Fab®	Isopropyl Alcohol	According to Manufacturer Spec	11-Aug-16
sBioMed, LLC	Steriplex SD Disinfectant Sporicide Wipes (Without Bleach)	Hydrogen Peroxide	According to Manufacturer Spec	20-May-15
Sterri-Matt®	Sterri-Matt® (Pure) Foam		Ready to Use	15-Dec-09
Spartan Chemical Company Inc.	Super HDQ Neutral®		256 to 1 dilution with water	11-Aug-16
PDI	Super Sani Cloth® Germicidal Disposable Wipes		Ready to Use	02-Sep-10
Reynard Health Supplies	Surface Sanitising Wipes		Ready to Use	18-Feb-09
3M	TB Quat Disinfectant	Ammonium Chloride	Ready to Use	09-Mar-18
Apollo Corporation	Therasol™		According to Manufacturer Spec	25-Sep-12
Medline Industries Inc.	Touch Screen Cleaning Wipes	70% Isopropyl Alcohol	Ready to Use	15-May-19
(Various)	Trix (Soap)		Ready to Use	27-Apr-10
Verna Care	Tuffie 5 Cleaning Wipes	Cocoalkyl Dimthylbenzyl Ammonium Chloride 0-0.5%	Ready to Use	19-Jul-12
Apollo Corporation	Turbo Clean™	Ethanol 9.5%	According to Manufacturer Spec	25-Sep-12
Clinell®	Universal Disinfectant Wipes		According to Manufacturer Spec	27-May-09



Approved Cleaners August 2019

PLEASE REFER TO IMPORTANT MESSAGES ON FINAL PAGE BEFORE APPLICATION OF ANY CLEANER/DISINFECTANT.

MANUFACTURER	PRODUCT	MAIN ACTIVE COMPONENT	CONCENTRATION	DATE TESTED
Clinell®	Universal Wipes		Ready to Use	12-May-19
True Blue Chemicals	Use All	Glycol 1-6%	According to Manufacturer Spec	28-Sep-12
Whiteley Medical	Viraclean®	Benzalkoniumchloride	According to Manufacturer Spec	19-Jul-12
Sealed Air Diversey Care	Virex® II 256		According to Manufacturer Spec	07-May-12
Sealed Air Diversey Care	Virex® Plus	Ammonium Chloride	According to Manufacturer Spec	19-Dec-18
Dupont™	Virkon®		1% Virkon (1:100 Water)	15-Dec-09
Wexford Labs	Wexcide 128		Minimum 128 to 1 dilution with water	11-Aug-16



Important Messages

CAUTION!

To avoid premature degradation of the surface coating never use an abrasive cloth, sponge or other similar device on the surface. Carefully follow manufacturer instructions for recommended exposure time of cleaner on surface. Wipe off cleaning product residue with a clean, soft cloth after cleaning and allow the surface to dry completely before use.

IMPORTANT INFORMATION:

Cleaners/disinfectants can cause discolouration/bleaching on some colours. To avoid discolouration, always check on a small area first.

Information contained here is purely related to the compatibility of cleaners with our materials and in no way relates to the efficacy of any cleaner for infection control or the safety of any cleaner for skin contact.

Testing info is specific to the products or formulations listed and in no way constitutes approval for other products from any brand listed without being tested.

All approvals are based on manufacturer recommended dilution with water unless otherwise specified and approval is void for any improper dilution.

If a particular cleaner/disinfectant is not listed, we advise to test the specific product in advance or to contact us to conduct testing.